

liveSCAN

Slide scanner and bright-field digital microscope



Your gateway to real-time
digital telepathology



LOGIBIOTECH

MAKING ROBOTICS. TAKING CARE.


Compact excellence for real-time digital telepathology

liveSCAN is an innovative **slide scanner** and **bright-field digital microscope** for **digital pathology**. Completely automated, fast, and reliable, **liveSCAN** enables the **scanning of slides** and the **visualization, analysis, archiving, retrieval, and sharing** of **digitized images**.



Our **liveSCAN** solution consists of a scanner-microscope unit with a loading capacity of 4 slides and 4 motorized objectives (4X, 10X, 20X, and 40X), an external PC with control software, and a 27-inch color LCD monitor.

Fast viewing and high quality images



LIVE MICROSCOPE

Immediate access to slide contents and “**live**” viewing with all objectives, even from **remote locations**, in real time, without the need for prior scanning.

DIGITAL MICROSCOPE


Virtual visualization and navigation of digitized slides with exceptional color fidelity and detail.

Predictive focus technology, high-performance objectives, and an advanced processing unit ensure the **highest quality images** and **fidelity of representation** at all times. The viewing experience is further optimized with **customizable features** that enhance image review and analysis.

Share, anytime and anywhere

for advanced collaborative diagnostics

With its features and sharing tools, **liveSCAN** supports **remote interaction** for **instant consultations**, allowing multiple users to **simultaneously navigate digital slides** in real time and **discuss cases** through **simple and intuitive collaboration options**.



A resource at the service of patients

The capability to share slides in real time with maximum security and diagnostic reliability, even remotely, represents **an additional tool to ensure the highest standards of patient care**.

Efficient workflow and optimal information management

liveSCAN integrates seamlessly with **laboratory information systems (LIS)**, creating a **collaborative work environment** within the digital pathology workflow, from case registration to slide management, archiving, and viewing of the final report.

The **software management system** allows for **efficient organization** of acquired slides, facilitating retrieval and utilization of information.



LOGIBIOTECH

LOGIBIOTECH srl

Registered Office

Via Fontana, 16 - 20122 Milano, Italy

Operating Headquarters

Via Emilia SS9 443 - 29010 Alseno (PC), Italy

www.logibiotech.com

e-mail: info@logibiotech.com

PEC: amministrazione@pec.logibiotech.com

VAT: IT10961670964

